

# KIMMON LASER SYSTEM

## IK3052R-BR


Wavelength : 325nm Power : 5mW

Transverse Mode : TEM multimode



Dimension (mm)

IK-B



## SPECIFICATIONS

Model	IK3052R-BR
Wavelength (nm)	325
Power (mW)	5
Transverse Mode	TEM multimode
Polarization	Random
Polarization Ratio	N/A
Noise, P-P @30kHz~2MHz (%)	< 8
Beam Diameter 1/e <sup>2</sup> (mm)	< 1.5* <sup>1</sup>
Beam Divergence (mrad)	< 0.8
Beam Pointing Stability ( $\mu$ rad) (At 25°C constant temperature)	$\leq \pm 25$
Warm Up Time (90% power) (minutes)	15
Power Stability (%) (At 25°C constant temperature)	$\leq \pm 2.0$ (4 hours)
Power Stability (10~40°C) (%)	< 20
Environmental Condition (Operation)	Temperature 10~40°C, Humidity $\leq 90\%RH^{*2}$
Environmental Condition (Storage)	Temperature -10~50°C, Humidity $\leq 90\%RH^{*2}$
Shock (with Kimmon packaging) (G)	20* <sup>3</sup>
Dimensions (W×H×L) (mm)	130×162×600
Weight (kg)	8.5
Laser Class	3B / IIIb

\*1 Measured at 100mm from output coupler

\*2 Non-condensing

\*3 Vertical direction

